

## 7F & 8F delivery system with hydrophilic coating

Catalog number	Nominal stent diameter (mm)	Nominal stent length (mm)	Use with introducer	Use with guidewire	Delivery system length (cm)
ZNSJ10x2035135	10	20	7F	0.035"	135
ZNSJ10x3035135	10	30	7F	0.035"	135
ZNSJ10x4035135	10	40	7F	0.035"	135
ZNSJ10x5035135	10	50	7F	0.035"	135
ZNSJ10x6035135	10	60	7F	0.035"	135
ZNSJ10x7035135	10	70	7F	0.035"	135
ZNSJ10x8035135	10	80	7F	0.035"	135
ZNSJ10x9035135	10	90	7F	0.035"	135
ZNSJ10x10035135	10	100	7F	0.035"	135
ZNSJ10x11035135	10	110	7F	0.035"	135
ZNSJ10x12035135	10	120	7F	0.035"	135
ZNSJ10x15035135	10	150	7F	0.035"	135

ZNSJ12x2035135	12	20	8F	0.035"	135
ZNSJ12x3035135	12	30	8F	0.035"	135
ZNSJ12x4035135	12	40	8F	0.035"	135
ZNSJ12x5035135	12	50	8F	0.035"	135
ZNSJ12x6035135	12	60	8F	0.035"	135
ZNSJ12x7035135	12	70	8F	0.035"	135
ZNSJ12x8035135	12	80	8F	0.035"	135

ZNSJ14x2035135	14	20	8F	0.035"	135
ZNSJ14x3035135	14	30	8F	0.035"	135
ZNSJ14x4035135	14	40	8F	0.035"	135
ZNSJ14x5035135	14	50	8F	0.035"	135
ZNSJ14x6035135	14	60	8F	0.035"	135
ZNSJ14x7035135	14	70	8F	0.035"	135
ZNSJ14x8035135	14	80	8F	0.035"	135

Caution:  
The standard length of delivery system is 135 cm, however 80 cm and 165 cm systems can be also ordered.  
To order, replace 135 by 80 or 165 in catalogue number – for example ZNSJ 7x303580, ZNSJ 7x3035165

Nominal stent diameters	4-14 mm
Nominal stent lengths	15-150 mm
Available delivery system length	80 cm, 135 cm, 165 cm
<b>NEW delivery system 6F</b>	For stents: diameter: 4-10 mm length: 15-150 mm
Delivery system 7F	For stents: diameter: 6-10 mm length: 15-150 mm
Delivery system 8F	For stents: diameter: 12-14 mm
Available guidewires	.035" or .018"
Delivery system type	OTW

# Jaquar

Self-expanding NiTi stent  
with delivery system

# Jaquar

Self-expanding NiTi stent  
with delivery system

**NEW delivery system!**

## Jaquar stent

- „One-wire“ technology
- Exceptional flexibility
- High radial force
- High resistant to:
  - kinking
  - crushing
  - fracturing

Caution:  
This publication could include technical or other inaccuracies or typographical errors.  
Changes are periodically added to the information on herein; these changes will be incorporated in the new editions of the publication. Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time. This publication is only a commercial information and cannot be handled the commercial offer in the meaning of the art.66 §1 of Polish Civil Code.

Self-expanding NiTi stent  
with delivery system

# Jaquar

**NEW delivery system!**

06/2011

**World's first  
„one-wire“ type NiTi stent  
since 2003!**

BALTON  
since  
1980

**BALTON**®

**BALTON**®

Vertrieb durch:  
**AVIDAL**  
DISTRIBUTION  
AVIDAL Distribution GmbH  
Argoallee 3, D-12527 Berlin  
Tel.: +49 30 67 89 69 84  
Fax: +49 30 67 89 69 85  
contact@avidal.net  
www.avidal.net/distribution

**BALTON**®

CE 1011

**AVIDAL**

